


Thalapump 20

Micropump™

Light and easy
to carry



Digital
programming
(screw driver free)



Clear display
of alarm
and speed



 **micrel**
Medical Devices

Designed
for patients comfort



Thalapump 20

Micropump™

The new MP Thalapump 20 benefits from the 20 years experience of Micrel in the field of Thalassaemia. This pump has been specially designed to fit the new standards of Iron Chelation for an outstanding patient care.

Specifications:

- Weight : 190 g with batteries (LR 03 AAAx6)
- Dimensions : 165x40x23 mm
- Battery autonomy : 70 full infusions minimum
- Accepted syringes: 5,10,20,30 ml
- Variable infusion speed: from 3 to 15 mm/h or 24 to 4 hours of infusion
- Prime facility
- Alarms and warnings:
 - End of travel
 - Occlusion
 - Low battery and depleted battery
 - Auto diagnostic

All information related to speed and alarms are clearly shown on the LCD display for a clear understanding of the infusion status.

Accessories:

The Thalapump™ 20 is provided with all the carrying accessories to provide patient comfort during the day.

CE
0120

Company ISO 9001/EN 46001
SGS U.K. certified
Federal law (U.S.) restricts
the MP Thalapump™ 20
for use by or on the order of a physician.



 **micrel**
Medical Devices

4 Ithakis Str., Pallini GR 15351 Athens Greece
Tel : +30 2106032333 & 4, Fax: +30 2106032335

www.micrelmed.com